



#### Welcome

Care Shop's premium brand, Curative, encompasses a range of products and equipment that aid the well being of residents.

- Pressure Area Care
- **Moving and Handling**
- **Specialist Seating**

The synonyms of 'Curative' are

#### Healing | Corrective | Therapeutic | Remedial

and that is exactly what our ranges do for your residents and all with maximum comfort and movement.

If you would like to know more about our ranges please call 0161 848 4000 or visit www.careshop.co.uk





#### Content

4	Bariatric Cushion
6	Static Castellated Cushion
8	Static High Risk Cushion
10	Comfort Static Mattress
12	Bariatric High Risk Mattress
14	Static High Risk Mattress
16	Static Memory Foam Mattress
18	Dynamic Pressure Care Mattress
20	Hybrid Pressure Care Mattress
22	Dynamic 8" Very High Risk Mattress

For your LOLER and PUWER inspections and servicing please call the DHS Service Team on 0161 848 4000

3













### Curative cs-NOR-C-444 **Bariatric Cushion**

- For High Risk residents
- Up to 50 stone (317kgs)
- Covers can be washed at 80 degrees
- Castellated foam surface
- Clinically proven

This specialist area of healthcare has additional levels of concerns to consider, which include a greater redistribution of the weight, additional support and airflow. To address these areas a new wave effect pattern has been developed to maximise comfort and to minimise the opportunity of pressure sores and ulcers.

The wave pattern changes across the length of the mattress to optimise pressure distribution in key areas. At both ends of the cushion they have a smaller depth of wave to enhance the density and comfort. As the cushion gets closer to the centre, the depth and width of the wave increase for greater redistribution of pressure and airflow.

This high-risk Bariatric cushion has been designed to be used with residents up to 50 stone (317kgs).

The multi-stretch material encapsulates the cushion using high-frequency welding to prevent fluid ingress and removes the opportunity of liquid contamination that would damage the cushion.

Its non-slip base will secure the cushion in place when residents stand or sit in a wheelchair or lounge chairs.

Fully compliant to British Standards and exceeds Crib 5 composite flame-retardant testing (BS7176), the Curative Bariatric cushion will provide optimum comfort and safety for your residents.

For your LOLER and PUWER inspections and servicing please call the DHS Service Team on 0161 848 4000

5













#### Curative cs-ECO-116-E Static Castellated Cushion

- For High Risk residents
- Castellated memory foam surface
- Fire retardant
- 2-Year warranty
- Clinically proven

Care Shops Curative static pressure care cushion has been developed to improve comfort whilst reducing pressure in the key areas of the gluteus maximus and lower back.

Utilising the latest innovations in temperature sensitive memory foams and high-density foam bases, our cushions will provide your residents with exceptional comfort.

With the castellation of the foam, the Curative Static Pressure Reduction Cushion is designed to flex and move with the residents. This dramatically reduces a constant single point of contact with the cushion, helping to minimise the risk of pressure sores and ulcers.

To further enhance the results, the cushion is encased with a multi-stretch material to re-distribute additional. pressure, reducing shear and friction sores.

Using high-frequency welding ensures all seams are firmly secured. This prevents fluid ingress and removes the opportunity of liquids contaminating and damaging the foams.

With its non-slip base, this static pressure cushion is perfect to use in wheelchairs and lounge chairs when required.

This cushion can be used for all levels of risk and completely confirms to all the British Standards, exceeding Crib 5 composite flame-retardant testing (BS7176). Healthcare environment needs.



#### **Technical Specification**

**Care Shop Product Code:** CS-VIS-115-3

Risk Factor: High Risk Fatigue Class:

**Maximum Patient Weight:** 17" (W) x 17" (L) x 3" (D) Dimension:

43cm (w) x 43cm (L) 7.5cm (D)

Less than 3 kgs **Cushion Weight:** 

210-230g (Options Available) **Cushion Cover Coating:** Exceeds minimum requirements **Fire Retardant Standards:** 

of Crib 5 Composite testing BS7176

Warranty: 2 Year











## Curative cs-vis-115-3 Static High Risk Cushion

- For High Risk residents
- 2-Year warranty
- Clinically proven
- Fire retardant

Care Shops Curative Static Pressure Care Cushion has the same great features as the CS-ECO-116-E 'Castellated' version but has been specifically designed for the High Risk residents.

It's made from the same materials that provide the great comfort and continues to use the high-end memory foam to reduce pressure in the key areas. As this is for High Risk residents the digitally engineered castellation is not required so the height is reduced to 3 inches (75mm)

The multi-stretch material encapsulates the cushion using high-frequency welding to prevent fluid ingress and removes the opportunity of liquid contamination that would damage the cushion.

Its non-slip base will secure the cushion in place when residents stand or sit in a wheelchair or lounge chairs.

With our 2-year warranty, we are confident that our Curative Static Pressure Cushion will provide the comfort and support your Healthcare environment needs whilst complying to British Standards and exceeding Crib 5 composite flame-retardant testing (BS7176).













### Curative ECOW/198/88/10 **Comfort Static Mattress**

- For Medium Risk residents
- Multi stretch material
- Fire retardant
- Crib 5 tested
- 2-Years warranty

The Curative Comfort Static Mattress has been developed for residents at medium risk of developing pressure sores.

This static mattress has been designed to evenly distribute the weight of a resident, providing comfort as well as facilitating pressure ulcer prevention.

The cover is made from a fire retardant, multi-stretch material that can be washed at high temperatures, reducing the risk of outbreaks.













# Curative <sub>CS-NOR-141</sub> Bariatric High Risk Mattress

- For High Risk residents
- Up to 50 stone (317kgs)
- Covers can be washed at 80 degrees
- Multi-stretch material
- None-slip base

The Curative Bariatric range has been developed to support the differing body shapes and weights from the heavier patient group. This specialist area of healthcare has additional levels of concerns to consider.

This includes a greater redistribution of the weight, additional support and airflow. To address these areas a new wave effect pattern has been developed to maximise comfort and to minimise the opportunity of pressure sores and ulcers.

The wave pattern changes across the length of the mattress to optimise pressure distribution in key areas. At both ends of the cushion they have a smaller depth of wave to enhance the density and comfort for the head and feet. As the mattress gets closer to the centre, the depth and width of the wave increase for greater re-distribution of pressure and airflow.

This high-risk Bariatric mattress has been designed to be used with residents up to 50 stone (317kgs). The multi-stretch material encapsulates the cushion using high-frequency welding to prevent fluid ingress and removes the opportunity of liquid contamination that would damage the cushion.

Its non-slip base will secure the mattress in place when residents stand or sit in a wheelchair or lounge chairs.

Fully compliant to British Standards and exceeds Crib 5 composite flame-retardant testing (BS7177), the Curative Bariatric cushion will provide optimum comfort and safety for your residents.













## Curative <sub>CS-BF-242</sub> Static High Risk Mattress

- For High Risk residents
- Castellated foam
- Multi-stretch material
- · Fits on a standard single bed
- Clinically proven

Our High Risk Curative mattresses has been designed to provide our customers with a high-quality low-cost solution for patients who need more support than the traditional mattress.

Utilising the technology created for our High and Very High Risk mattresses, we have used a similar design to maximise comfort and the redistribution of pressure.

The high-density combustion modified, ether foam (CME) maintains exceptional support for the resident whilst providing excellent airflow around the mattress to moderate the temperature for additional comfort.

As with all our mattresses, they are encased with a multistretch, waterproof material to enhance the redistribution of pressure and reduce the potential of sores and ulcers. The mattress covers are hand sown in UK with a zip flap cover to reduce the opportunity of fluid ingress.

Our static mattresses exceed all the BS standards and comply with even the strictest of infection prevention control protocols.













# Curative CS-VIS-M-112 Static Memory Foam Mattress

- For High Risk residents
- Castellated memory foam
- Multi-stretch material
- Fire retardant
- 5-Year warranty

Our Static Pressure Care Mattresses are created using the latest Combustion Modified Ether (CME) foam technology to provide maximum stability and breathability.

The Castellated foam surface creates multiple contact points to maximise the pressure redistribution. It also helps with the airflow around the mattress assisting the control of temperature and improving the comfort and ability to sleep.

The specially designed U core surrounds the Castellated foam and provides additional rigidity and support.

At the head and foot ends of the bed, the castellated foam is shaped smaller to support a protect the most vulnerable areas, where the larger castellation in the middle, maximises the support of the body.

As with all our mattresses, they are encased with a multistretch, waterproof material to enhance the redistribution of pressure and reduce the potential of sores and ulcers. Using high-frequency welding, the foam is encapsulated with this Crib 5 composite tested material and prevents fluid ingress and a non-slip base is added as standard.

Our static mattresses exceed all the BS standards and comply with even the strictest of infection prevention control protocols and come with a 5-year Guarantee.













# Curative cs-5F-121 Dynamic Pressure Care Mattress

- For High to Very High Risk residents
- Digital control panel
- Extremely quiet motor
- 2-Year warranty
- CPR function

When residents develop bed pressure sores and ulcers, standard mattresses will hinder the healing process and often inflict greater damage to the skin.

Our Dynamic mattresses can provide 24 hours of rapid pressure relief to make your high to very high-risk patients as comfortable as possible whilst their sores and ulcers heal.

With one of the quietest air pumps on the market today, your patients will be able to sleep, whilst the dynamic mattress goes through its rapid inflation and deflation cycles. This not only maximises the recuperation of their symptoms but provides a comfortable environment.

With a simple interface on the control panel, it can be quickly used to adjust pressure settings to match individual patients' weights and comfort requirements.

Each of the 22 Thermoplastic Polyurethane (TPU) cells is rigorously tested to ensure the bed can be used by people of weights up to 28 stone, in any environment.

The rapid deflate functionality is an important feature for when CPR is required to be performed. The pull cord on the mattress will deflate the cells in seconds leaving a flat surface against the high-density foam to perform this emergency procedure.













## Curative CS-HYB-MO (MATTRESS) HYB-121/PUMP (PUMP) Hybrid Pressure Care Mattress

- For High to Very High Risk residents
- Combined dynamic and static mattress
- CPR function
- Ultra-soft heel area
- Digital control panel
- 2-Year warranty

Utilising the combined innovations of our dynamic and castellation mattresses we have developed the Curative Hybrid Mattress. It enhances the comfort and recovery time of very high-risk residents with up to Grade 4 pressure ulcers.

To reduce heal damage, our patient pending capillary wave design redistributes pressure directly from the heal to increase the surface contact area to include the calf. This

reduces shear and friction forces whilst improves blood flow to improve the overall comfort of the resident.

The head and body area of the mattress uses our clinically proven dynamic solution that maintains pelvic stability whilst helping the carers in the moving and handling procedures.

In the event of an emergency response, the Curative hybrid mattress provides a good surface to perform CPR. Simply detach the mattress from the pump and the foam filled cells (inflated or deflated) will provide the necessary rigidity to execute the chest compressions.

We are confident that the Curative Hybrid Mattress will be a perfect addition to your home especially with its 2-year warranty across the mattress and pump and its universal use.

#### **Curative Hybrid Digital Pump** (Sold Separately)

Our digital pump is one of the quietest on the market and has the ability to inflate in 5 kg increments. This enables this mattress to be used for light weight frail individuals through to residents up to 39 stone.













# Curative cs-DC-118 Dynamic 8" Very High Risk Mattress

- For Very High Risk residents
- CPR function
- Digital control panel
- 2-Year warranty

The Curative Dynamic 8" Mattress is one of Care Shop's Premium Mattresses and designed for residents who are at a VERY HIGH RISK to bed sores or who have existing tissue damage (pressure ulcer grade 1 to 4).

This dynamic mattress has been developed to adjustment to each resident's individual body shape and position, for optimum pressure care.

This mattress has been created to provide the resident with maximum comfort. This is achieved by the alternating action of each cell that provides lower interface pressures and gentle capillary stimulation.

Our digital pump is one of the quietest on the market and has the ability to inflate in 5 kg increments. This enables this mattress to be used for light weight frail individuals through to residents up to 28 stone.



Dynamic | Castellated | Hybrid | Static | Bariatric



