



# Welcome

Care Shop's premium brand, Curative, encompasses a range of products and equipment that aid the well being of residents.

- **Pressure Area Care**
- **Moving and Handling**
- **Specialist Seating**

The synonyms of 'Curative' are

**Healing | Corrective | Therapeutic | Remedial**

and that is exactly what our ranges do for your residents and all with maximum comfort and movement.

If you would like to know more about our ranges please call 0161 484 4000 or visit [www.careshop.co.uk](http://www.careshop.co.uk)



## Content


4	About Curative
6	Universal Sling
8	Classic Universal Sling
10	Toilet Sling
12	Hammock Spacer Sling
14	Stand aid/Back support belt




Sonic welding for additional strength



Laser guided stitching



BSEN 10535



Wide range of Materials



# Curative Moving and Handling range by Care Shop

- **Laser guided stitching of high quality materials to maximise strength and durability.**
- **Available in 4 different sizes and 5 different materials**
- **Suitable for high temperature washing up-to 90 degrees**
- **High Definition woven label that will last as long as the sling**
- **Loop and Clip fixings across the range**
- **1 Year guarantee across the range**
- **Exceeds the British Standard EN 10535**

Care Shop's moving and handling range, Curative, has been designed to maximise comfort and safety for both the resident and health care professional.


Our entire range is manufactured in the UK and surpasses the strict British Standards guidelines.

To maximise your investment, Delivered Health Solutions service and Installation teams will support you and your equipment.


They specialise in on-site LOLER, PUWER and PAT testing across all moving, handling and pressure area care equipment.

Training is also available on all moving and handling equipment from Care Shop .


For your LOLER and PUWER inspections and servicing please call the DHS Service Team on 0161 848 4015




Through-feed dignity hoop







Optional head sections






Wide range of Materials

Colour	Fabric
Quilted Teal Polyester	

Size	Weight Limit	Product Code (Loop)	Product code (Clip)	Binding
Small	< 220kg / 35 st	CUR/UNI/POLY/S	CUR/UNI/POLY/S/CL	
Medium	< 220kg / 35 st	CUR/UNI/POLY/M	CUR/UNI/POLY/M/CL	
Large	< 220kg / 35 st	CUR/UNI/POLY/L	CUR/UNI/POLY/L/CL	
X-Large	< 220kg / 35 st	CUR/UNI/POLY/XL	CUR/UNI/POLY/XL/CL	

Please contact Customer Services for additional product codes, which will vary according to the specifications required.  
Bespoke / Made to measure service available.

Other colours available	Fabric
Grey Mesh Polyester	
Teal Solid Polyester	
Black Spacer Fabric	
Blue Parasilk	



BS EN 10535



Flame Retardant



Warranty




High Temp Wash



Made in the UK

# Curative Universal Sling



- **General purpose sling**
- **Easy to apply in a chair, bed or from floor positions.**
- **Available with or without head section**
- **Through-feed dignity hoop**
- **High Definition woven label that will last as long as the sling**

The Universal Sling is one of our most popular products in the Curative range, due to its versatility, comfort and design.


The appropriately cut leg sections, combined with our soft binding, provide a comfortable lift for both the resident and the health care provider.

The cut of the sling and versatile sizing provides pelvic stability throughout the lifting procedure.

Optional head sections are available to ensure patients requiring support for the neck and head regions are accommodated.

A wide range of materials are available, including quick-drying mesh fabric for bathing and water therapy processes.

For your LOLER and PUWER inspections and servicing please call the DHS Service Team on 0161 848 4015



Laser guided stitching,



HD Colour Woven label



Coloured coded loop

Colour	Fabric
Quilted Teal Polyester	

Size	Weight Limit	Product Code (Loop)	Product code (Clip)	Binding
Small	< 220kg / 35st	CUR/CLA/UNI/POLY/S	N/A	
Medium	< 220kg / 35st	CUR/CLA/UNI/POLY/M	N/A	
Large	< 220kg / 35st	CUR/CLA/UNI/POLY/L	N/A	
X-Large	< 220kg / 35st	CUR/CLA/UNI/POLY/XL	N/A	

Please contact Customer Services for additional product codes, which will vary according to the specifications required.  
Bespoke / Made to measure service available.

Other colours available	Fabric
Grey Mesh Polyester	
Teal Solid Polyester	
Black Spacer Fabric	
Blue Parasilk	



BS EN 10535



Flame Retardant



Warranty




High Temp Wash



Made in the UK



# Curative Classic Universal Sling



- **Wider than our Universal sling with additional support**
- **Separate leg areas**
- **Can be applied underneath both legs in a hammock position.**
- **Easy to apply in a chair, bed or from floor positions.**

Our Classic Universal sling has been designed with separate leg areas, providing simpler functionality.

The easy fit sling offers great positioning, supporting the entire body during the lifting process, especially around the hip area.

It can be applied from either a lying or sitting position and is suitable for some amputee patients, subject to the necessary assessment.

Our laser guided stitching enhances the strength of our slings, surpassing the EN10535 manufacturing standard.

The high quality, anti-shrink material can be washed at high temperature, ensuring an aesthetically pleasing product throughout the life of the sling.

For your LOLER and PUWER inspections and servicing please call the DHS Service Team on 0161 848 4015






Appropriately cut leg sections,







Coloured coded loop







HD Colour Woven label

Colour	Fabric
Quilted Teal Polyester	

Size	Weight Limit	Product Code (Loop)	Product code (Clip)	Binding
Small	< 220kg / 35 st	CUR/TOILET/POLY/S	CUR/TOILET/POLY/S/CL	
Medium	< 220kg / 35 st	CUR/TOILET/POLY/M	CUR/TOILET/POLY/M/CL	
Large	< 220kg / 35 st	CUR/TOILET/POLY/L	CUR/TOILET/POLY/L/CL	
X-Large	< 220kg / 35 st	CUR/TOILET/POLY/XL	CUR/TOILET/POLY/XL/CL	

Please contact Customer Services for additional product codes, which will vary according to the specifications required.  
Bespoke / Made to measure service available.

Other colours available	Fabric
Grey Mesh Polyester	
Teal Solid Polyester	
Black Spacer Fabric	
Blue Parasilk	



BS EN 10535



Flame Retardant



Warranty




High Temp Wash



Made in the UK

# Curative Toilet Sling



- **Through-feed hoops help with the residents dignity when transferring them to the toilet**
- **Ultra comfortable arm padding creating a relaxing lift, without unwanted compression**
- **Residents must have some upper body strength to use this sling**



Our premium toilet sling has been designed to ensure the toilet transfer process is made easier, and clothing removal is greatly facilitated.

The Curative toilet sling can be used with most loop spreader bar hoists and is available with head supports, if required.

The through-feed dignity hoop, ensures ease of application, without the need for additional webbing.

We use a premium Velcro® webbing belt, that provides a strong adhesion, with maximum comfort. An additional belt and buckle can be added as an upgrade to provide greater support for the resident.

For your LOLER and PUWER inspections and servicing please call the DHS Service Team on 0161 848 4015



Colour	Fabric
Black Spacer Fabric	

Size	Weight Limit	Product Code (Loop)	Product code (Clip)	Binding
Small	< 220kg / 35 st	CUR/HAMM/SPC/S	CUR/HAMM/SPC/S/CL	
Medium	< 220kg / 35 st	CUR/HAMM/SPC/M	CUR/HAMM/SPC/M/CL	
Large	< 220kg / 35 st	CUR/HAMM/SPC/L	CUR/HAMM/SPC/L/CL	
X-Large	< 220kg / 35 st	CUR/HAMM/SPC/XL	CUR/HAMM/SPC/XL/CL	

Please contact Customer Services for additional product codes, which will vary according to the specifications required.  
Bespoke / Made to measure service available.

Other colours available	Fabric
Grey Mesh Polyester	
Teal Solid Polyester	
Quilted Teal Polyester	
Blue Parasilk	



BS EN 10535



Flame Retardant



Warranty




High Temp Wash



Made in the UK

# Curative Hammock Sling



- **Designed for residents with severe immobility/amputees and potential risk of falling**
- **Apply from the bed or floor**
- **Colour coded webbing to help accuracy when lifting**
- **Unique soft strap protector**


The Curative Hammock sling is manufactured in various materials and sizes to accommodate the differing needs of the residents and health care professionals.

The slings can be washed at high temperatures and when manufactured in the Spacer material offers additional pressure reducing qualities.

The Hammock sling is available with loops or clips and can be used with most ceiling track products and mobile hoists.

The soft strap protector, which is unique to our sling designs, provides security to the patient in the most vulnerable of areas.

For your LOLER and PUWER inspections and servicing please call the DHS Service Team on 0161 848 4015



Colour	Fabric
Quilted Teal Polyester	

Size	Weight Limit	Product Code (Loop)	Product code (Clip)	Binding
Small	< 220kg / 35 st	CUR/SA/POLY/S	N/A	
Medium	< 220kg / 35 st	CUR/SA/POLY/M	N/A	
Large	< 220kg / 35 st	CUR/SA/POLY/L	N/A	
X-Large	< 220kg / 35 st	CUR/SA/POLY/XL	N/A	

Please contact Customer Services for additional product codes, which will vary according to the specifications required.  
Bespoke / Made to measure service available.

Other colours available	Fabric
Grey Mesh Polyester	
Teal Solid Polyester	
Black Spacer Fabric	
Blue Parasilk	



BS EN 10535



Flame Retardant



Warranty




High Temp Wash



Made in the UK





# Curative Stand aid/Back support belt

- **Aids your residents to stand and sit whilst supporting their spine**
- **Strengthening stays create additional rigidity and support to the lower back area**
- **Padded for additional comfort**
- **To be used with residents that have a high level of body tone**

The Curative Stand Aid/Back Support belt has been developed in our British manufacturing facilities, using high-quality materials and innovative technology.

Our design provides complete support for the patient, and peace of mind for the health care professional

The slings strengthening stays provide rigidity in the back region, maximising the support in the lower and central spinal areas.

We have added additional padding to provide extra comfort to the resident and the soft strap protector provides security to the patient in the most vulnerable of areas.

For your LOLER and PUWER inspections and servicing please call the DHS Service Team on 0161 848 4015



Universal | Toilet | Hammock | Stand Aid | Bespoke  
LOLER and PUWER Servicing



Care Shop | Unit B Circle Court | Trafford Park | Warren Bruce Road | Manchester, M17 1LB

[www.careshop.co.uk](http://www.careshop.co.uk) | 0161 484 4000 | [sales@deliveredhealthsolutions.co.uk](mailto:sales@deliveredhealthsolutions.co.uk)