



## The Team

Our nationwide service and Installation teams specialise in onsite LOLER, PUWER and PAT testing across all moving, handling and pressure area care equipment.

All our employees are Dementia Friends have been specifically trained by manufacturers to install, test and repair their equipment. This helps to realise a strong return on investment for each asset, increasing the inlife use by up to 57%.

We adhere to the strict guidelines set out by the CQC, HSE and MHRA to ensure your organisation complies with the required standards and safety to protect you, your team and the residents.



## **Asset Management**

Slings and Accessories | Through Floor & platform lifts | Ceiling Track | Bath Lifters Bath Hoists Standing & Raising | Aids Hoists Profiling Beds | Wheelchairs | Pressure Area Care Mattresses | Shower Trolleys

The ability to track assets throughout your organisation provides complete visibility of their location, condition and quantity. Our team of engineers will visit your home to identify each asset with their unique serial number. The asset is then digitally catalogued detailing their specific room location and condition.

If the asset needs a repair, service or inspection, we will detail the required action and supply this information to your management team. When agreed, we will perform the necessary work to bring the asset back to a compliant standard and record the details.

We can also create a calendar reminder to alert you when your assets need to be re tested to comply with the LOLER and PUWER standards as part of the Asset Management service.



## LOLER and PUWER Servicing & Repair Installation and Training | Asset Management

Care Shop | Unit B Circle Court | Trafford Park | Warren Bruce Road | Manchester, M17 1LB careshop.co.uk

0161 848 4000

sales@deliveredhealthsolutions.co.uk